

an analyti**chem** company





# The Portable Microbial Air Sampler for Cleanroom Monitoring

The SAS range has been at the fore front of environmental monitoring for over 40 years. The distinctive yellow samplers are used by many of the world's leading pharmaceutical companies to demonstrate GMP compliance within their manufacturing facilities.

Designed and developed in Italy, the SAS Super range of portable microbial air samplers offer the user an easy to use, robust sampler for viable particle monitoring. Typical applications include monitoring within controlled spaces, such as clean rooms, operating theatres or laboratories, as well as office and outdoor environments.

The battery powered units will sample approximately 40,000 litres on a full charge and are easy to wipe down and clean. Each sampler is supplied complete with carry case, charger and calibration certificate.

## **Applications**

The SAS portable sampler can be used for many environmental monitoring applications, the key uses are:

- · Routine clean room monitoring
- Grade A monitoring of isolator cabinets
- Operating theatres

- Offices and HVAC systems
- · External composting sites
- Investigative work

#### **Features & Benefits**

#### Multi-point impaction sampling

Developed from multi stage sieve samplers, the SAS utilises the proven multi-point impaction method to collect airborne particles. 219 holes of 1mm diameter direct air flow on to the agar surface of a standard plate.

#### **Battery powered**

The nickel metal hydride battery pack will provide approximately 40,000 litres worth of sampling on a single charge. The rapid charger will have

the unit fully recharged in approximately 4 hours, allowing a significant number of samples to be completed per day.

#### Stainless steel sampling head

The 316 stainless steel Drilled Head can be steam autoclaved between uses or decontaminated using proprietary disinfectants or alcohol. Additional heads can be purchased so ensuring a ready to go head is always available.

#### **Delay start**

The SAS Super has a selectable delay start of 1 - 60 minutes at one minute specificity up to 10 minutes, followed by increments of 5 minutes up to 60 minutes. This allows the user to program the sampler and exit the room being monitored before the sample is taken. A true at rest sample is taken avoiding any impact from the operator.

#### Multi-mode sampling

If a prolonged sample to mimic a process is required, the SAS Super's Multi mode facility is ideal. A standard 1 cubic metre sample can be divided into smaller samples with a time delay between each.

#### Last sample memorised

The SAS Super makes it easy to repeat the same sample as it memorises the last sample set. Simply switch on and press Start to repeat the same sample volume.

#### **USB**

The SAS air sampler has been modernised with up-to-date USB technology to increase ease of use and user convenience:

- Easier data management via pdf download
- Data is downloaded by USB-C memory stick

### **Technical Data**

Product Specifications		
Sampling Rate:	SAS Super 100 – 100 litres/minute SAS Super 180 - 180 litres/minute	
Sample Size:	1 up to a maximum of 2000 litres operator selectable – 8 presets plus 8 user definable	
One m³ run time:	100L version = 10 mins, 180L version = 5 mins, 33 secs	
Construction Materials:	Powder coated polyurethane resin & stainless steel (full details available)	
Consumables Used:	Contact plates or 90mm Petri dishes – must be defined at ordering	
Power Supply:	Rechargeable Nickel Metal Hydride (without memory effect)	
Battery Capacity:	2450 mAh (approximately than 40,000 litres of air sampling)	
Battery Recharge Time:	Charge for at least 4 hours	
Sampler Dimensions (mm):	290 x 105 x 110 (H x W x D	
Sampler Weight (kg):	With stainless steel drilled head – 2.2kg	
Carry Case Dimensions (mm):	455 x 368 x 166 (H x W x D)	
Total Weight (kg) (inc sampler, case, charger):	With stainless steel drilled head – 5.3kg	

Ordering Information		
Part Number	Description	Packaging
3.1129	SAS Super 100 USB air sampler, Contact, s/s drilled head	Per Item
3.1132	SAS Super 100 USB air sampler, Petri, s/s drilled head	Per Item
3.1130	SAS Super 180 USB air sampler, Contact, s/s drilled head	Per Item
3.1133	SAS Super 180 USB air sampler, Petri, s/s drilled head	Per Item
Optional Extra	as	
3.1070	SAS Tripod 1.44m floor standing	Each
3.2005	SAS Stand – stainless steel	Each
3.1069	SAS Tripod Telesco 2.9m Floor standing	Each
2.1252	IQ/OQ Manual	Each
3.1065	SAS Daily Heads - Contact	Case of 40
3.1066	SAS Daily Heads – 90mm Petri	Case of 40
3.1092	SAS Drilled Head S/S Contact	Each
3.1093	SAS Drilled Head S/S Petri	Each

#### **Terms & Conditions**

All SAS air samplers are supplied as complete kits, ready to use, including calibration certificate and a 24 months return to base warranty. We reserve the right to amend specifications without notice. E&OE



7 & 8 Launton Business Centre, Murdock Road, Bicester, OX26 4XB, United Kingdom