

# CERTIFICATE OF CALIBRATION



**Certificate Number:** 04/600

**Company:** ABC PHARMA LTD  
ABC HOUSE  
MAIN ROAD  
NEW TOWN  
XXI IXX

**Date Certified:** 4 JULY 2006

**Model:** SAS SUPER 180                      **Serial Number:** NN-NN-NNNN

**Status:** RECALIBRATION                      **Instrument Used:** TA5 A

**Report Ref:** SR xxxxx

## CALIBRATION RESULTS

**\*Nominal Air Flow Rate:**                      **180 litres/minute**  
(This confirms the unit is sampling  
within  $\pm 3\%$  of the nominal air flow)

**Time Cycle (200 litres):**                      **66 seconds**

**Calibrated By:** xxxxxxxx                      **Signed:**

**We recommend that this instrument is re-calibrated before: 4 July 2007**

**\* This is to certify that the above instrument has been calibrated in accordance with our methodology. See reverse for details.**

**Form: Calib12 Issue 2 10/02**

