

# QUALITY CONTROL CERTIFICATE

## REDIPOR IRRADIATED PREPARED MEDIA



**Product Code :** 104130  
**Description :** Tryptone Soya Agar  
 90mm Petri Dishes  
  
**Batch Number :** 220612  
**Date of Manufacture :** 22 June 2006  
**Expiry Date (see Note) :** 22 December 2006  
**pH (T=7.1 ±0.2) :** 7.1  
**Irradiation Reference :** R20606 0569      **Mean Irradiation Dose (T=20-35 kGy) :** 29.6 kGy

*Positive Control organisms to show a ratio of 0.7 and above (Inoculum level 10-100 cfu)*

<i>Organism</i>	<i>Reference</i>	<i>Productivity Ratio</i>
Staphylococcus aureus	ATCC 6538P	1.0
Staphylococcus epidermidis	ATCC 14990	0.8
Escherichia coli	ATCC 8739	0.9
Pseudomonas aeruginosa	ATCC 9027	0.9
Bacillus subtilis	ATCC 6633	1.7
Candida albicans	ATCC 10231	1.0

*Positive Control organism to show growth*

Aspergillus niger	ATCC 16404	Growth & Black Spores
-------------------	------------	-----------------------

*Comments:*

*Signed By:*

*Authorised Signatory*

*Date:* 3 July 2006

**NOTE:** Expiry Date only valid if packs stored unopened at Ambient Room Temperature not exceeding 25°C