

# Newsletter

Autumn/Winter 2008



www.cherwell-labs.co.uk

## Redipor Barrier Pack - the advantages are clear!

Redipor Barrier Pack provides the best solution for contact plates and settle plates in gassed isolators.

The inner layer of wrap is an easy to tear, heat sealed, transparent film, impermeable to H<sub>2</sub>O<sub>2</sub> vapour. Outside the seal it has sufficient film to make tearing with gauntlets easier and a hole for hanging. It is supplied as an irradiated triple wrap with the outer layers of transparent polythene for robust protection until it is required in the isolator.

**redipor**  
BARRIER PACK



Any of the irradiated Redipor plates can be supplied in the Barrier Pack presentation on request.

Please contact our Sales Department ([sales@cherwell-labs.co.uk](mailto:sales@cherwell-labs.co.uk)) for more information.

## Redipor - batch numbers & product codes

Cherwell Laboratories has embarked on an exciting journey to modernise our IT systems with the aim of providing a more efficient service to our customers. We will also be updating our batch numbering system and introducing a new product code system across the whole business.

From **January 2009, Redipor batch numbers will contain a year code.** We will write to customers individually with the details as well as putting them on our website: [www.cherwell-labs.co.uk](http://www.cherwell-labs.co.uk)

**During 2009 our product codes will change to an eight digit number** to allow for future growth. We will operate both old and new codes in parallel for 12 months to allow all customers to order with confidence while we update their information.

We will inform all customers before we implement this change and hope to publish the new codes in our April 2009 price list.

Please contact us if you have any specific questions about these changes.

## Christmas!

You should hopefully have received within the last week details about our deliveries over the Christmas holidays. If not a copy of the information is on our website at:

[www.cherwell-labs.co.uk/news-10-08.php](http://www.cherwell-labs.co.uk/news-10-08.php)

This year, as in previous years, we will be making donations directly to charities rather than sending Christmas cards.


**May we take this opportunity to wish seasons greetings to all our customers!**

## More SAS accessories!

The new SAS stainless steel stand is designed to improve the stability of an SAS sampler when placed vertically for sampling.

If your samplers are prone to being knocked over and damaged the SAS Stand could be the perfect answer. It is easy to fit and even allows the operator to program the unit whilst fitted.

For further information, please contact Sales.



SAS

November  
2008



**Cherwell Laboratories Ltd**

7 & 8 Launton Business Centre Murdock Road BICESTER OX26 4XB  
T: 01869 355500 F: 01869 355545 E: [sales@cherwell-labs.co.uk](mailto:sales@cherwell-labs.co.uk)