Certificate of Analysis REDIPOR PREPARED MEDIA



Part Number 4.0005 Old code 100014

Product Description Blood Agar Base + 5% Sheep blood 90mm Petri dish

single vent clear 18ml fill. Flow wrap. Sleeve of 10.

Batch reference 25030742
Work Order Number W2507709
Date of Manufacture 07 March 2025
Expiry Date 11 April 2025

pH (T=7.1 to 7.5) 7.5 Sterility Test Pass (Absence of growth following 72 hours at 30 - 35°C)

Positive Control organisms to show growth (Inoculum level 10-100 cfu)

Organism	Reference	Result
Bacillus spizizenii	ATCC 6633	Growth
Escherichia coli	ATCC 8739	Growth
Pseudomonas paraeruginosa	ATCC 9027	Growth
Staphylococcus aureus	ATCC 6538	Growth
Candida albicans	ATCC 10231	Growth
Clostridium sporogenes	ATCC 19404	Growth
Staphylococcus epidermidis	ATCC 12228	Growth

Released By: Ela Szokalo Date: 14 March 2025

All our prepared media products are manufactured and tested at our site in Bicester, Oxfordshire under the control of our ISO9001 : 2015 approved quality management system. The generation of this certificate confirms all sterilisation and performance criteria have been achieved

NOTE: Expiry Date only valid if packs stored unopened at Ambient Room Temperature not exceeding 25°C Cherwell Laboratories Ltd, 7 & 8 Launton Business Centre, Murdock Road, Bicester, 0X26 4XB