Redipor® Triple Wrap Pack


Cherwell Triple-Wrapped Irradiated Plates Packaging

The diagram below presents Cherwell's Redipor® Triple wrap pack packaging, comprising of three layers with two layers of flow wrap and a large sealed polybag. This is a representation of our standard triple wrap that can be applied across our range of plated product.

The illustration demonstrates the sterility assurance with each layer, following exposure to gamma radiation. It allows the plates to be easily handled without the risk of introducing contaminants, while transporting product through each stage of a facility.

Maintaining conditions for the agar and providing consistent, reliable prepared media that delivers confidence in sampling.





Cardboard Box

- Sealed on all edges
- Irradiated
- Labelled with Irradiation Indicator and Product Information

Inner Plastic Layer

- Easy Tear Incision on closure
- Labelled with Irradiation Indicator and Product Information
- Material: PET/LLDPE

Outer Plastic Layer

- Easy Tear Incision on closure
- Material: PET/LLDPE

Outer Polythene Bag

- Sealed
- Material: LDPE