

Growth Media Type - Sabouraud Dextrose Agar (SDA)

Organisms & morphology: Sabouraud Dextrose Agar will support growth of xx. A typical set of control organisms and the phenotypic identification associated with these are:

Fungal	Morphology		
	Edge	Shape	Colour
Candida albicans	Smooth	Regular	Cream
Aspergillus brasiliensis	Irregular	Irregular	White to black



Email:
sales@cherwell-labs.co.uk



Website:
www.cherwell-labs.co.uk



Telephone:
01869 355500

Growth Media Type - Sabouraud Dextrose Agar (SDA)

Regulatory references: SDA is a Harmonised Pharmacopoeia medium and as such is mentioned in the BP/EP/JP, as well as the US Pharmacopoeia.

EP Chapter Ref 2.6.12 Microbial examination of non-sterile products: Microbial enumeration tests

Recommended Culture Media	Property	Micro-organisms	Incubation Time & Temp
Sabouraud Dextrose Agar	Growth promotion	Candida albicans Aspergillus brasiliensis	20-25°C for < 5 days



Email:
sales@cherwell-labs.co.uk



Website:
www.cherwell-labs.co.uk



Telephone:
01869 355500